

# Adjustable Bed - (Portable) 4 Function


PRODUCT CODE

## ABF-P4F



- This bed is fully portable, folding onto a frame utilising its 4 lockable casters for ease of movement.
- The bed can be stored in this configuration, thereby saving considerable space over conventional beds.
- Assembly or dismantling takes about 10 minutes.

It has 4 functions by hand held controller:



1. Raise & lower bed



2. Tilt back rest



3. Raise & lower knee break



4. Tilt bed ends (Trendelenburg/reverse trendelenburg function)

- Width: 90cm
- Length: 218cm
- Height Adjustment: 36 - 80cm
- Weight: 60kg
- Max User Weight: 180kg



Distributed By

**Mobility Rentals & Sales**  
**1300 460 070**

Unit 22 / 170 North Road,  
 Woodridge 4114, QLD Australia  
[www.mobilityrentals.com.au](http://www.mobilityrentals.com.au)