

PRODUCT DATA SHEET Adult Multi-Functional Wheelchair	Created Date 30/01/2011	Version 1.0
	Author RDS	
PRODUCT CODE WAMAP-46S	File Name WAMAP-46S.pdf	
	File Location	...MRSA\CompanyGraphics(NEW)\DataSheets-MRSA



- ADJUSTABLE HEIGHT ARM RESTS
- COMFORT PADDED SEAT & BACK REST
- ADJUSTABLE BACK REST
- ELEVATING LEG RESTS
- ADJUSTABLE SEAT
- ALUMINIUM FRAME
- RECLINING
- TILT-IN-SPACE



Multi functional wheelchair, designed for people with Cerebral Palsy

Typical users

Tilt-in-space wheelchairs are often used by adults who have:

- Progressive medical conditions, moderate to severe motor involvement.
- Poor head and trunk control, with limited ability to sit upright.
- Limitations in range of motion, especially hips, or with curvature of the spine (scoliosis or kyphosis).
- Pain or pressure while sitting, requiring contoured seat to relieve pressure.
- Weakness, lack of endurance, fluctuating muscle tone, seizure or sleep disorders.

91cm	100cm	49cm	26cm	40cm	15cm	40cm	46 cm	52cm	68 - 74 cm	16 > 8.6kg	135kg

Distributed By



Unit 22 / 170 North Road,
Woodridge 4114, QLD Australia
www.mobilityrentals.com.au